OXYSTAT

COLORIMETRIC ASSAY FOR THE QUANTITATIVE DETERMINATION OF PEROXIDES
IN EDTA PLASMA, SERUM AND OTHER BIOLOGICAL FLUIDS
CAT. NO. BI-5007 12 x 8 TESTS

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FOR RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC PROCEDURES

rev.no. 170504 (replacing 150804)

This kit was developed and manufactured by:
Biomedica Medizinprodukte GmbH & Co KG, A-1210 Wien, Divischgasse 4
Tel. +43/1/291 07 45, Fax +43/1/291 07 85, E-mail export@bmgrp.com

CONTENT / INHALT / SOMMAIRE / CONTENIDO

- 1) ENGLISH..... 3**
- 2) DEUTSCH ... 7**
- 3) FRANCAIS...11**
- 4) ESPANIOL... 15**

*Additional Information about our products is available on our homepage
Zusätzliche Information zu unseren Produkten ist auf unserer Homepage erhältlich
Des informations complémentaires sur notre gamme de produits sont disponibles sur notre site Internet
Encontrará información adicional sobre nuestros productos en nuestra página web*

www.bmgrp.com

1) INTRODUCTION

Cells and tissues are sensitive to oxidative stress, caused by the formation of free radicals. If not deactivated by antioxidants, organic peroxides and hydroperoxides are the first reaction products between cellular constituents and free radicals or other reactive oxygen derivatives.

The determination of the oxidative status / oxidative stress is essential in today's medical research and diagnostics. Methods used so far were either expensive (HPLC), or detected only degradation products of polyunsaturated fatty acids, like TBARS (thiobarbituric acid reactive substances).

The Biomedica OxyStat assay measures the total peroxide concentration of a sample, utilizing a quick and simple assay procedure. Results show a direct correlation between free radicals and circulating biological peroxides and thus allow the characterization of the oxidative status in biological samples.

Possible Indications

- Cardiovascular disease
- Atherosclerosis
- Sepsis
- Carcinogenesis
- Inflammatory processes
- Neurodegenerative processes

2) CONTENTS OF THE KIT

CONT	KIT COMPONENTS	QUANTITY
PLATE	Microtiter strips in strip holder	12 x 8 tests
SOLA	Solution A, Sample buffer, ready to use	1 x 25 ml
SOLB	Solution B, Reaction buffer, ready to use	1 x 1 ml
SOLC	Solution C, Enzyme solution, ready to use	1 x 50 µl
SOLD	Solution D, Reconstitution solution, ready to use	1 x 5 ml
CTRL 1+ 2	Controls, lyophilised. The concentration after reconstitution is stated on the label	2 vials lyophilised
CAL	Calibrators, lyophilised. The concentration after reconstitution is stated on the label	3 vials lyophilised
STOP	STOP solution, sulphuric acid, ready to use	1 x 6 ml

3) ADDITIONAL MATERIAL ADDED TO THE KIT

- Instruction manual for use
- Protocol sheet

4) EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Precision pipettes calibrated to deliver 5 µl – 5000 µl and disposable tips
- ELISA reader for absorbance at 450 nm
- Incubator for 37°C

5) REAGENTS AND SAMPLE PREPARATION

We recommend EDTA-plasma as sample type because in serum a time dependent increase in peroxide concentration is observed. When serum is chosen, please make sure that during preparation of serum a time period no longer than 30 min at room temperature is allowed for clotting. Store EDTA or serum samples at -20°C.

Heparinized plasma should not be used for this assay, because it often precipitates upon storage and thus changes the results.

Heparin-plasma, lipemic or haemolytic samples may give erroneous results.

Cloudy samples should be centrifuged at least 5 minutes at 5000xg before use in the assay. All samples should be mixed well before assaying.

Reconstitute as follows:

- Dissolve one CAL (calibrator) and each CTRL (control) in 250 µl solution D and leave at room temperature (18-26°C) for five minutes, mix well. Reconstituted CAL and CTRL can be stored aliquoted and frozen at -20°C up to one month.
- The following volume of ABC-reaction-mix is sufficient for 40 tests (wells) and should be used as guidance. However, the volume should be adjusted to the respective number of samples. Prepare the ABC-reaction-mix immediately before the assay

by mixing:
5 ml solution A + 100 µl solution B + 5 µl solution C.

6) PRINCIPLE OF THE ASSAY

The peroxide concentration is determined by reaction of the biological peroxides with peroxidase and a subsequent colour-reaction using TMB as substrate. After addition of a stop solution, the coloured liquid is measured photometrically at 450 nm.

A calibrator is used to calculate the concentration of circulating biological peroxides in the sample (one-point calibration).

7) ASSAY PROTOCOL

All reagents and samples must be at room temperature (18-26°C) before use in the assay.
Mark position for CTRL/CAL/SAMPLE (Control/Calibrator/Sample) on the protocol sheet.
Take microtiter strips out of the plastic bag and mark as appropriate. Store unused strips in the bag.
Add 10 µl CTRL/CAL/SAMPLE (Control/Calibrator/Sample) into respective well.
Add 100 µl SOLA (Solution A Sample buffer) into each well.
Measurement 1: Determine absorption with an ELISA reader at 450 nm.
Add 100 µl of the freshly prepared ABC-reaction-mix into all wells (see section 5) Reagents and Sample preparation).
Incubate 15 minutes at 37°C.
Add 50 µl STOP (Stop solution) into each well.
Measurement 2: Determine absorption with an ELISA reader at 450 nm.

Assay procedure for automated pipetting systems

- Adjust volumes of the reaction solutions and samples according to the analyser used.
- If the analyser is incapable of adding the stop-solution, perform kinetic measurements at 405 nm or at 580- 620 nm without stopping the reaction.
- Measurement time: 15 minutes at 20- 37°C.

8) CALCULATION OF RESULTS

The differences between measurement 1 and 2 are proportional to the peroxide concentrations of the samples.

- For each CTRL/CAL/SAMPLE (Control/Calibrator/Sample) subtract the OD-values of measurement 1 from the OD-values of measurement 2 (= Δ OD)
- A single-point calibration is done using the CAL (Calibrator). The OD-value is proportional to its concentration, which is stated on the label. This concentration is stated as H₂O₂-equivalents (µmol/l).
- The concentrations of controls and samples are calculated according to following formula:

$$[\mu\text{mol} / \text{l}] \text{ sample} = \frac{\Delta \text{OD sample} \times [\mu\text{mol} / \text{l}] \text{ calibrator}}{\Delta \text{OD calibrator}}$$

9) ASSAY CHARACTERISTICS

Normal range:	Reference values from apparently healthy persons with no documented disease and medication: EDTA-plasma: <400 µmol/l Serum: <350 µmol/l Median: 372 µmol/l Each laboratory should establish own normal values.
Linearity:	Up to 600 µmol/l l
Sample volume:	10 µl human EDTA plasma, serum and other biological fluids
Incubation time:	15 min
Detection Limit:	7 µmol/l

10) PRECISION

Intra-assay (n = 12)	
Mean ($\mu\text{mol/l}$)	221
SD ($\mu\text{mol/l}$)	6.9
CV (%)	3.1

Inter-assay (n = 12)	
Mean ($\mu\text{mol/l}$)	221
SD ($\mu\text{mol/l}$)	11.3
CV (%)	5.1

11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date.
- Protect reagents from direct sunlight.
- Avoid foaming when mixing reagents.
- Stop solution should be added to the plate in the same order as the reaction-mix solution.

12) PRECAUTIONS

All test components of human source were tested against HIV-Ab, HCV-Ab and HBsAg; and were found negative. Nevertheless, they should be handled and disposed as if they were infectious.

Liquid reagents contain $\leq 0.1\%$ Proclin 300 as preservative.

Avoid contact with skin and mucous membrane. Proclin 300 is not toxic in concentrations used in this kit. It may cause allergic skin reactions – avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke or apply cosmetics where reagents are used.
- Avoid all contact with reagents by using gloves.
- Sulfuric acid is irritating to eyes and skin. Flush with water if contact occurs. Avoid contact with skin and mucous. Irritations are possible. – Flush with water after contact!!

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BI-5007 OXYSTAT

ASSAY PROTOCOL AND CHECKLIST

- Bring all reagents to room temperature (18-26°C).
- Prepare reagents and samples as instructed.
- Bring unused and prepared components to the storage temperature mentioned in the package insert.
- Take microtiter strips out of the bag and mark positions on the protocol sheet.

- 1. Add 10 µl CTRL/CAL/SAMPLE (control/calibrator/sample) into each well, except blank.
- 2. Add 100 µl SOLA (solution A sample buffer) into each well.
- 3. Measurement 1: Determine absorption with an ELISA reader at 450 nm.**
- 4. Add 100 µl freshly prepared ABC reaction mix into each well.
- 5. Incubate for 15 minute at 37°C.**
- 6. Add 50 µl STOP (Stop solution) into each well.
- 7. Measurement 2: Determine absorption with an ELISA reader at 450 nm.**