


## NT-proBNP ELISA (Cat.No. SK-1204)

### For the Determination of NT-proBNP (1-76) in Human Samples

#### **ASSAY CHARACTERISTICS**

<b>Method</b>	Sandwich ELISA, HRP/TMB, 12x8-well strips
<b>Antibodies/Standard</b>	Capture antibody: polyclonal anti human NT-proBNP antibody, specifically binding to amino acids 31-57 of proBNP. Detection antibody: polyclonal anti human NT-proBNP antibody, conjugated to peroxidase, specifically binding to amino acids 8-29 of proBNP. Standard material: synthetic human NT-proBNP (1-76).
<b>Sample type</b>	Human serum, EDTA plasma
<b>Standard range</b>	0-640 pmol/l (0, 10, 40, 160, 640 pmol/l)
<b>Conversion factor</b>	1 pmol/l = 8.475 pg/ml (MW: 8.475 kDa) <i>refers to NT-proBNP (1-76) that is detected by the ELISA</i>
<b>Sample volume</b>	50 µl / well – direct measurement
<b>Sensitivity</b>	LOD (0pmol/l + 3SD): 3.0 pmol/l; LLOQ: 3.3 pmol/l
<b>Specificity</b>	This assay recognizes endogenous (natural) and recombinant human NT-proBNP (1-76).
<b>Incubation</b>	3 h /30 min – room temperature

#### **Typical standard curve of Biomedica NT-proBNP ELISA:**



**NT-proBNP ELISA, #SK-1204 – ASSAY & PERFORMANCE CHARACTERISTICS**
**Serum values from apparently healthy donors:**

<b>Serum (n=70)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>5.8</b>
Mean	9.8
Maximum	63.9
Minimum	1.1
Percentil 95%	32.5
Percentil 5%	2.0

**EDTA plasma values from apparently healthy donors:**

<b>EDTA plasma (n=28)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>5.6</b>
Mean	9.5
Maximum	32.0
Minimum	0.6
Percentil 95%	27.9
Percentil 5%	0.8

**Serum values from unselected hospital panel:**

<b>Serum (n=117)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>46.5</b>
Mean	84.5
Maximum	621.1
Minimum	0.0
Percentil 95%	22.1
Percentil 5%	0.8

**EDTA plasma values from unselected hospital panel:**

<b>EDTA plasma (n=40)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>28.6</b>
Mean	57.9
Maximum	489.2
Minimum	0.3
Percentil 95%	157.7
Percentil 5%	1.2

**Serum values of a heart failure panel NYHA classification II-IV measured with the Biomedica NT-proBNP ELISA (#SK-1204)**

<b>Serum (n=66)</b>	<b>NT-proBNP (pmol/l)</b>
mean	47.5
max	399.3
min	1.1

**NT-proBNP ELISA, #SK-1204 – ASSAY & PERFORMANCE CHARACTERISTICS**
**NYHA classification subgroups:**

measured with the Biomedica NT-proBNP ELISA (#SK-1204)

Data on serum values of a heart failure panel **NYHA classification II**


Serum (n=43)	NT-proBNP (pmol/l)
mean	24.8
max	75.6
min	1.1

Data on serum values of a heart failure panel **NYHA classification III**

Serum (n=16)	NT-proBNP (pmol/l)
mean	56.2
max	264.3
min	4.7

Data on serum values of a heart failure panel **NYHA classification IV**

Serum (n=7)	NT-proBNP (pmol/l)
mean	166.8
max	399.3
min	84.1

**MATRIX COMPARISON**
**Correlation of serum and plasma samples from apparently healthy individuals**


## **PERFORMANCE CHARACTERISTICS**

### **Spike Recovery in human serum samples:**

**The mean recovery of recombinant NT-proBNP in serum samples (n=4) when spiked with two different concentrations of human recombinant NT-proBNP is 99 and 108%.**

#### Experiment:

Recovery of spiked samples was tested by adding 2 concentrations of human recombinant NT-proBNP (1-76) to 4 different human serum samples.

Data showing spike/recovery of human serum samples:

<b>Matrix</b>		<b>Serum</b>			<b>S/R (%)</b>	
		<b>0</b>	<b>80</b>	<b>320</b>	<b>80</b>	<b>320</b>
<b>Sample ID</b>	#1	1.6	88.5	347.2	109	108
	#2	0.1	74.2	360.2	93	113
	#3	1.7	80.9	338.7	99	105
	#4	4.4	81.6	346.5	96	107
				<b>Mean (%)</b>	<b>99</b>	<b>108</b>

### **Spike Recovery in human EDTA plasma samples:**

**The mean recovery of recombinant NT-proBNP in EDTA plasma samples (n=4) when spiked with two different concentrations of human recombinant NT-proBNP is 93 and 94%.**

#### Experiment:

Recovery of spiked samples was tested by adding 2 concentrations of human recombinant NT-proBNP (1-76) to 4 different human EDTA plasma samples.

Data showing spike/recovery of human EDTA plasma samples:

<b>Matrix</b>		<b>EDTA plasma</b>			<b>S/R (%)</b>	
		<b>0</b>	<b>80</b>	<b>320</b>	<b>80</b>	<b>320</b>
<b>Sample ID</b>	#1	0.8	73.4	276.8	92	87
	#2	7.7	72.8	279.0	91	87
	#3	9.3	77.2	324.1	97	101
	#4	19.1	74.9	309.5	94	97
				<b>Mean (%)</b>	<b>94</b>	<b>93</b>

### **Dilution Linearity/Parallelism:**

#### **Serum (n=3):**

**mean 117%. Range: 112-119%**

(samples spiked with recombinant NT-proBNP)

**mean 82%. Range: 79-86%** (endogenous NT-proBNP)

#### **EDTA plasma (n=2):**

**mean 84%. Range: 83-84%**

(samples spiked with recombinant NT-proBNP)

**mean 84%. Range: 79-80%** (endogenous NT-proBNP)

Experiment 1 - Dilution Linearity in human serum samples:

Dilution linearity was assessed by spiking serum samples with recombinant NT-proBNP (1-76) (320 pmol/l) and by diluting the samples with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit):

Data showing the dilution of recombinant NT-proBNP (1-76) in serum samples:

<b>Serum Sample ID</b>	<b>Reference NT-proBNP (pmol/l)</b>	<b>Spiked NT-proBNP (pmol/l)</b>	<b>Dilution 1+1</b>	
			<b>NT-proBNP (pmol/l)</b>	<b>R (%)</b>
#1	1.7	338.7	202.2	119
#2	0.1	360.2	202.4	112
#3	8.0	278.5	165.7	119
			<b>Mean (%)</b>	<b>117</b>

Experiment 2 - Parallelism in human serum samples:

Dilution linearity (parallelism) was assessed by diluting serum samples containing endogenous NT-proBNP with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit).

Data showing the dilution of endogenous NT-proBNP (1-76) in serum samples:

<b>Serum Sample ID</b>	<b>Reference NT-proBNP (pmol/l)</b>	<b>Dilution 1+1</b>	
		<b>NT-proBNP (pmol/l)</b>	<b>R (%)</b>
#1	44.2	17.4	79
#2	51.3	20.4	80
#3	61.7	26.5	86
			<b>Mean (%)</b>
			<b>82</b>

**EDTA plasma (n=2):** **83-84%** (spiked with recombinant NT-proBNP)  
**80-85%** (endogenous NT-proBNP)

Experiment 1 - Dilution Linearity in human plasma samples:

Dilution linearity was assessed by spiking EDTA plasma samples 1+1 with recombinant NT-proBNP (1-76) (640 pmol/l) and by diluting the samples with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit):

Data showing the dilution of recombinant NT-proBNP (1-76) in EDTA plasma samples:

<b>EDTA plasma Sample ID</b>	<b>Reference NT-proBNP (pmol/l)</b>	<b>Spiked NT-proBNP (pmol/l)</b>	<b>Dilution 1+1</b>	
			<b>NT-proBNP (pmol/l)</b>	<b>R (%)</b>
#1	36	347	145	84
#2	0.8	277	115	83
			<b>Mean (%)</b>	<b>84</b>

Experiment 2 – Parallelism in human plasma samples:

Dilution linearity was assessed by diluting EDTA plasma samples containing endogenous NT-proBNP with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit).

Data showing the dilution of endogenous NT-proBNP (1-76) in EDTA plasma samples:

EDTA plasma Sample ID	Reference	Dilution 1+1	
	NT-proBNP (pmol/l)	NT-proBNP (pmol/l)	R (%)
#1	44.2	17.4	79
#2	51.3	20.4	80
<b>Mean (%)</b>		<b>84</b>	

Intra-assay (With-In Run) & Inter-assay precision (In-Between)

**Intra-assay (n=3) ≤ 4%, Inter-assay (n=8) ≤ 7%**

Experiment:

Intra-assay: 2 samples of known concentrations were tested 3 times in 1 assay by 1 operator.

Inter-assay: 2 samples of known concentrations were tested 8 times in 2 assays by different operators.

Data showing intra-assay and inter-assay precision:

Intra-assay (n=3)	Sample 1	Sample 2	Inter-assay (n=8)	Sample 1	Sample 2
Mean (pmol/l)	60.2	35.2	Mean (pmol/l)	52.1	108.1
SD (pmol/l)	2.0	0.9	SD (pmol/l)	1.7	7.9
CV (%)	4	3	CV (%)	3	7

The limit of detection (LOD):

The LOD is defined as the mean value of the back calculated concentration plus 3 times the standard deviation. The LOD of the NT-proBNP ELISA is 3 pmol/l.

The lower limit of quantification (LLOQ):

The lower limit of quantification is defined as the accuracy of the back calculated concentrations and shall not exceed ±25% (acc. to ICH (Ref. 1)).

For the NT-proBNP ELISA the LLOQ was determined as 3.3 pmol/l.

Specificity:

This assay recognizes endogenous (natural) and recombinant human NT-proBNP (1-76).

Calibration:

This immunoassay is calibrated against recombinant human NT-proBNP (1-76) peptide.

## SAMPLE CHARACTERISTICS

### Sample stability:

We recommend separating EDTA plasma or serum by centrifugation as soon as possible, e.g. 20 min at 2,000 x g, preferably at 4°C (2-8°C). EDTA-plasma or serum can be stored at 4°C (2-8°C) up to two days. For long term storage, aliquot the acquired EDTA plasma or serum samples and store them at -25°C or lower.

Samples can be subjected to 3 freeze-thaw cycles.


Serum and EDTA samples can be stored for at least 2 years at -80°C.

### Freeze/thaw of samples containing endogenous NT-proBNP

The mean recovery of serum and EDTA plasma sample concentrations stressed by 3 F/T cycles is 94%. The mean CV of sample concentrations (not stressed and stressed up to 3 times by freeze-thaw cycles) is 4%.

NT-proBNP concentrations of samples after 3 freeze-thaw cycles:

<b>Sample ID</b>	<b>NT-proBNP [pmol/l]</b>					<b>CV (%)</b>	<b>R (%) 3 F/T vs ref</b>
	<b># of F/T cycles</b>	<b>Ref</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>Mean</b>	
Serum1	70	67	67	69	68	1	98
Serum2	15	13	15	14	14	5	98
Serum3	215	185	191	188	194	6	87
Serum4	141	129	129	126	131	4	90
Serum5	49	47	48	47	48	2	96
Plasma1	203	191	199	187	195	3	92
Plasma2	282	257	271	248	264	5	88
Plasma3	68	64	79	69	70	8	101
<b>Mean (%)</b>						<b>4</b>	<b>94</b>



**Validation**

The assay is fully validated according to ICH Q2 (R1), Ref. (1).

**References**

(1) CPMP/ICH/381/95 - ICH Topic Q2 (R1) „Validation of Analytical Procedures: Text and Methodology“ including:

ICH Q2A “Text on Validation of Analytical Procedures”

ICH Q2B “Validation of Analytical Procedures: Methodology”

**Available on our website**

[www.bmgrp.com](http://www.bmgrp.com)


[Enzyme immunoassay for the quantitative determination of NT-proBNP \(1-76\) in human serum \(cat.no. SK-1204\)](#)

Protocol booklet, Validation Data, MSDS, Citations

## NT-proBNP ASSAY COMPARISON

Comparison between NT-proBNP Assays: Biomedica #SK-1204 versus Roche Cobas



### 1. Comparison of human serum samples (n=66): NYHA (II-IV)



Comparison of serum values of a heart failure panel **NYHA classification II-IV** measured with the Biomedica NT-proBNP ELISA and the Roche Cobas assay (*values shown both in pmol/l and in pg/ml*)

<b>n=66</b>	<b>Biomedica</b> <i>pmol/l</i>	<b>Roche</b> <i>pmol/l</i>	<b>Biomedica</b> <i>pg/ml</i>	<b>Roche</b> <i>pg/ml</i>
<b>mean</b>	47.5	33.9	402.5	287.3
<b>max</b>	399.3	398.7	3384.0	3378.9
<b>min</b>	1.1	1.1	9.3	9.3


## 2. Comparison of NYHA (IV) samples (n=7)



Comparison of serum values of a heart failure panel **NYHA classification IV** measured with the Biomedica NT-proBNP ELISA and the Roche Cobas assay (*values shown both in pmol/l and in pg/ml*)

n=7	Biomedica pmol/l	Roche pmol/l	Biomedica pg/ml	Roche pg/ml
<b>mean</b>	166.8	130.9	1413.6	1109.2
<b>max</b>	399.3	398.7	3384.0	3379.0
<b>min</b>	84.1	28.1	712.7	238.4


### 3. Comparison of NYHA (III) samples (n=16)



Comparison of serum values of a heart failure panel **NYHA classification III** measured with the Biomedica NT-proBNP ELISA and the Roche Cobas assay (*values shown both in pmol/l and in pg/ml*)

<b>n=16</b>	<b>Biomedica</b> <i>pmol/l</i>	<b>Roche</b> <i>pmol/l</i>	<b>Biomedica</b> <i>pg/ml</i>	<b>Roche</b> <i>pg/ml</i>
<b>mean</b>	56.2	44.5	476.2	377.9
<b>max</b>	264.3	200.1	2239.9	1696.0
<b>min</b>	4.7	4.3	39.8	36.6

#### 4. Comparison of NYHA (II) samples (n=43)



Comparison of serum values of a heart failure panel **NYHA classification II** measured with the Biomedica NT-proBNP ELISA and the Roche Cobas assay (values shown both in pmol/l and in pg/ml)

n=43	Biomedica pmol/l	Roche pmol/l	Biomedica pg/ml	Roche pg/ml
<b>mean</b>	24.8	14.0	210.1	119,3
<b>max</b>	75.6	44.7	640.7	379,1
<b>min</b>	1.1	1.1	9.3	9,5