Novel ELISA for the quantification of human leucine-rich α-2 glycoprotein (LRG) in serum and plasma

Elisabeth Gadermaier¹, Jacqueline Wallwitz¹, Gabriela Berg^{, 2}, Gottfried Himmler¹

¹The Antibody Lab GmbH, Vienna, Austria ²Biomedica Medizinprodukte GmbH, Vienna, Austria

SUMMARY AND CONCLUSION

Special features of the LRG sandwich ELISA

- High precision
- Characterized antibodies
- Wide assay range
- Includes: 7 standards, 2 controls
- Correlates well with other assays
- Full validation package
- For serum and plasma

The

Lab

Antibody

MEDICA

Good sample stability

INTRODUCTION

Leucine-rich α -2 glycoprotein 1 (LRG) is a highly glycosylated 38 kDa protein. It belongs to the leucine-rich repeat protein family whose members are involved *e.g.* in proteinprotein interactions, cell signaling, or cell adhesion.

It is described as proangiogenic factor and increased levels are reported in various diseases where neovascularinvolved. zation is For instance, it was shown that TNF-α LRG induced by promotes angiogenesis that then contributes to aberrant formation in the bone subchondral bone of osteoarthritis patients. A recent study highlights the clinical inflamusefulness of the matory marker LRG for the evaluation of disease activity arthritis. rheumatoid IN Furthermore, LRG may also a biomarker for serve as rheumatoid arthritis disease activity during IL-6 blockade treatment.

LEUCINE-RICH ALPHA-2-GLYCOPROTEIN (LRG) VALUES IN DISEASE PANELS



			LRG [µg/ml]	
Samples	n	Mean	Range	Median
Apparently healthy controls	18	27.7	19.2 - 40.2	27.5
Rheuma cohort	18	53.5	33.0 - 79.0	51.0
Arthrose cohort	16	39	12.0 - 89.0	33.0
Kidney disease	16	54.7	36.0 - 90.0	47.0

LEUCINE-RICH ALPHA-2-GLYCOPROTEIN (LRG) ELISA CHARACTERISTICS

Principle of the detection of human LRG



The LRG ELISA is a sandwich-based immunoassay. Antibodies were characterized by epitope mapping of linear epitopes with microarray technology and by the determination of binding kinetics with biolayer interferometry. Plates are coated with a peptide-specific polyclonal sheep anti-human LRG antibody binding to a linear epitope in the N-terminal region of LRG. Multiple linear epitopes recognized by the polyclonal detection antibody are distributed over the whole LRG sequence and are located in the N- and C-terminus, as well as within the leucinerich repeats. Both antibodies bind to LRG with low dissociation

METHODS

Here, we show the development, characterization and validation of a new LRG ELISA (cat. no. BI-LRG, Biomedica, Vienna, Austria).

The assay is optimized and validated for human serum and plasma samples according to international quality guidelines. Urine and cell-culture supernatants are compatible with this ELISA.

STABILITY OF LEUCINE-RICH ALPHA-2-GLYCOPROTEIN (LRG)

Benchtop stability of human samples

Bench top stability of serum samples



Benchtop stability of human serum samples containing endogenous LRG. Samples were stored for 3 hours at room temperature (rt) or overnight (o.n.) at 4°C and compared with the reference resulting in less than 4% measurement alterations.

Freeze-thaw stability of human samples

Freeze thaw stability of samples

25,0

rate constants.

LRG assay protocol

100 µl sample (prediluted) / standard / control			
(5µl sample volume required)			
2h @ 18-26°C			
5x washing			
100µl detection antibody			
2h @ 18-26°C			
5x washing			
100 µl substrate			
30 min @ 18-26°C (dark)			
50 µl stop solution read OD (450 nm)			

Samples, standards (0 - 64 ng/ml) and controls are pipetted into the antibody coated plate. After a washing step the labelled detection antibody is incubated for 2 hours followed by an additional washing step. Thereafter, TMB substrate is added for 30 minutes, followed by stopping the reaction and OD measurement at 450 nm.

LRG ELISA validated for serum and plasma



Assay performance, sample matrix comparison as well as sample measurements of apparently healthy and diseased human subjects were performed.



Freeze-thaw stability (F/T) (x-axis) of three human samples containing endogenous serum LRG. Samples are stable for at least 5x freeze-thaw cycles.

LITERATURE

- MicroRNA-497 elevation or LRG1 knockdown promotes osteoblast proliferation and collagen synthesis in osteoporosis via TGF-β1/Smads signalling pathway.
 Gu Z et al., J Cell Mol Med. 2020; 24(21):12619-12632.
- Leucine-rich alpha2-glycoprotein as the acute-phase reactant to detect systemic juvenile idiopathic arthritis disease activity during anti-interleukin-6 blockade therapy: A case series. Shimizu M et al., Mod Rheumatol. 2017; 27:833–837.
- Usefulness of Serum Leucine-Rich Alpha-2 Glycoprotein as a Disease Activity Biomarker in Patients with Rheumatoid Arthritis. Ha, You Jung et al., jkms. 2014; 29(9):1199–1204.

LRG was measured in serum and plasma samples prepared from 10 apparently healthy donors. A mean CV values of 10% was observed between all measured sample matrixes which lies within the recommended validation quality guidelines. The assay can be used for all indicated samples matrices.

Gabriela Berg:

CONTACT

gabriela.berg@bmgrp.com