# FGF23 (INTACT)

(EN) ELISA FOR THE QUANTITATIVE DETERMINATION OF HUMAN  
INTACT FGF23 IN EDTA PLASMA, HEPARIN PLASMA, AND CITRATE PLASMA  
Cat. No. BI-20700 . 12 x 8 TESTS

(DE) ELISA ZUR QUANTITATIVEN BESTIMMUNG VON HUMANEM,  
INTAKTEM FGF23 IN EDTA PLASMA, HEPARIN PLASMA UND CITRAT PLASMA  
Kat. Nr. BI-20700 . 12 x 8 TESTE

FOR RESEARCH USE ONLY  
NOT FOR USE IN DIAGNOSTIC PROCEDURES

This ELISA is optimized and validated for human plasma samples. Serum and other sample types are compatible with this ELISA (for more information please visit our website).

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Not for Sale in the United States!

Rev.no. 220524 (replacing 190612)

This kit was developed and manufactured by:  
Biomedica Medizinprodukte GmbH, A-1210 Wien, Divischgasse 4  
Tel. +43/1/29107 45, Fax +43/1/29107 6389, e-mail [info@bmgrp.com](mailto:info@bmgrp.com)

## **CONTENT / INHALT**

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Detailed information on the FGF23 (intact) human ELISA, e.g. assay validation data, sample matrix comparisons, and stability data is available on our website.

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***[www.bmgrp.com](http://www.bmgrp.com)***

## 1) INTRODUCTION

FGF23 (fibroblast growth factor 23) is a member of the fibroblast growth factor family and controls phosphate and vitamin D homeostasis. The full-length protein comprises 251 amino acids including a 24 amino acid signal peptide. The N-terminal FGF region of FGF23 is separated from the C-terminal region by a proteolytic cleavage site. A fraction of FGF23 is proteolytically processed between arginine179 and serine180 to generate N-terminal and C-terminal fragments. Therefore, the main forms of FGF23 present in human circulation are hormonally intact FGF23 and inactive N-terminal and C-terminal fragments. FGF23 binds to FGF receptor 1c (FGFR1c) with its N-terminal region, while the C-terminal region interacts with the co-receptor  $\alpha$ -Klotho to confer high-affinity binding to the receptor. FGFR1c and  $\alpha$ -Klotho are expressed in the distal nephron and the parathyroid gland. Co-receptor independent signaling of FGF23 has been described for other FGFRs, which are expressed in a variety of tissues. The main source of FGF23 are osteocytes in the bone.

Areas of interest:

### Bone Diseases

Autosomal dominant hypophosphatemic rickets (ADHR), autosomal recessive hypophosphatemic rickets (ARHR), X-linked hypophosphatemic rickets (XLH), McCune-Albright syndrome/fibrous dysplasia, hypophosphatemic rickets and hyperparathyroidism (HRHPT), drug-induced rickets/osteomalacia, hyperphosphatemia, hyperostosis-hyperphosphatemia syndrome, fibrous dysplasia, neurofibromatosis, osteoglophonic dysplasia (OGD), linear nevus sebaceous syndrome (LNSS)/epidermal nevus syndrome (ENS), Jansen-type metaphyseal chondrodysplasia

### Cancer

Tumor-induced osteomalacia, familial tumor calcinosis

### Cardiovascular Disease

Heart failure, left-ventricular hypertrophy, vascular calcification

### Endocrine disorders

Primary hyperparathyroidism, secondary hyperparathyroidism

### Kidney Disease

Chronic kidney disease (CKD), chronic kidney disease-mineral and bone disorder (CKD-MBD), chronic kidney disease and cardiovascular disease, acute kidney injury

### Metabolic disorders

Dyslipidemia, type 2 diabetes, anemia

## 2) CONTENTS OF THE KIT

CONTENTS	KIT COMPONENTS	QUANTITY
PLATE	Recombinant monoclonal anti-human FGF23 antibody precoated microtiter strips in strip holder, packed in aluminium bag with desiccant	12 x 8 tests
WASHBUF	Wash buffer concentrate, 20x, natural cap	1 x 50 ml
ASYBUF	Assay buffer, red cap, ready to use	1 x 10 ml
STD	Standards (0; 50; 100; 200; 400; 800; 1600 pg/ml), recombinant human FGF23 in human plasma, white caps, lyophilized	7 vials
CTRL	Controls A and B, yellow caps, lyophilized, exact concentration see labels	2 vials
CONJ	Conjugate (monoclonal mouse anti-human FGF23-HRP), amber cap, ready to use	1 x 7 ml
SUB	Substrate (TMB solution), blue cap, ready to use	1 x 13 ml
STOP	Stop solution, white cap, ready to use	1 x 7 ml

## 3) ADDITIONAL MATERIAL IN THE KIT

- 1 self-adhesive plastic film
- Quality control protocol
- Protocol sheet
- Instruction for use

#### 4) MATERIAL AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Precision pipettes calibrated to deliver 10 µl, 50 µl, 100 µl, 400 µl, and disposable tips
- Distilled or deionised water
- Plate washer is recommended for washing, alternative multichannel pipette or manifold dispenser
- Refrigerator with 4°C (2-8°C)
- Microplate reader capable of measuring absorbance at 450 nm (with correction wavelength at 630 nm)
- EP-tubes for sample dilution if pre-dilution is necessary
- Graph paper or software for calculation of results

#### 5) REAGENTS AND SAMPLE PREPARATION

All reagents as supplied in the kit are stable at 4°C (2-8°C) until the expiry date stated on the label of each reagent.

##### Sample preparation/dilution:

Collect venous blood samples by using standardized blood collection tubes. Perform plasma or serum separation by centrifugation according to supplier's instructions of the blood collection devices. **Assay the acquired samples immediately or aliquot and store at -25°C or lower.** Lipemic or haemolyzed samples may give erroneous results. Samples are stable for at least four freeze-thaw cycles. Samples should be mixed well before assaying. We recommend duplicates for all values.

**Samples with values above Standard 7 (STD7) can be diluted 1:11 (1+10), e.g. 10 µl sample + 100 µl ASYBUF, and tested again.**

The kit includes sufficient assay buffer (ASYBUF) for a 1:11 dilution of 40 samples (in duplicates). Additional ASYBUF can be ordered using cat.no. BI-20700-ASYBUF.

Due to the instability of the intact FGF23 molecule samples should be collected, processed, aliquoted and stored at -25°C or below as quickly as possible.

For further information on sample stability please visit our website [www.bmgrp.com](http://www.bmgrp.com) (see Validation Data) or contact our customer service by e-mail [info@bmgrp.com](mailto:info@bmgrp.com) or by phone +43/ 1/ 29107-45.

##### Reconstitution/handling:

**WASHBUF (wash buffer):** Dilute the concentrate 1:20, e.g. 50 ml WASHBUF + 950 ml distilled or deionized water. Crystals in the buffer concentrate will dissolve at room temperature. The undiluted WASHBUF is stable at 4°C (2-8°C) until the expiry date stated on label. The diluted WASHBUF is stable up to one month at 4°C (2-8°C). Only use diluted WASHBUF when performing the assay.


**STD (standards) + CTRL (controls):** Pipette 400 µl of distilled or deionized water into each vial. Leave at room temperature (18-24°C) for 15 min. Vortex gently. The exact concentration is printed on the label. Reconstituted STD and CTRL can be stored at -25°C or lower until expiry date stated on the label. STDs and CTRLs are stable for four freeze-thaw cycles.

#### 6) PRINCIPLE OF THE ASSAY

The FGF23 (intact) human ELISA kit is a sandwich enzyme immunoassay that has been optimized and validated for the quantitative determination of intact FGF23 in human plasma samples. Serum and other sample types are compatible with this ELISA. For further information on serum measurements please see page 6/7 or visit our website (see Validation Data).

In a first step, standard/control/sample and conjugate (monoclonal mouse anti-human FGF23-HRP) are pipetted into the wells of the microtiter strips, which are pre-coated with a recombinant monoclonal anti-human FGF23 antibody. FGF23 present in the standard/control/sample binds to the pre-coated antibody in the well and forms a sandwich with the conjugated anti-human FGF23-HRP antibody. In a washing step all non-specific unbound material is removed. In a second step, the substrate (TMB, tetramethylbenzidine) is pipetted into the wells. The enzyme-catalyzed color change of the substrate is directly proportional to the amount of intact FGF23 present in the sample. This color change is detectable with a standard microplate reader. A dose response curve of the absorbance (optical density, OD at 450 nm) using the values obtained from the standards versus the standard

concentration is generated. The concentration of intact FGF23 in the sample is determined directly from the dose response curve.

## 7) ASSAY PROTOCOL

All reagents and samples must be at room temperature (18-24°C) before use in the assay.

Mark position for STD/CTRL/SAMPLE (standard/control/sample) on the protocol sheet.

Take microtiter strips out of the aluminium bag. Store unused strips with desiccant at 4°C (2-8°C) in the aluminium bag. Strips are stable until expiry date stated on the label.


1. Pipette 50 µl STD/CTRL/SAMPLE (standard/control/sample) in duplicate into respective well.
2. Add 50 µl CONJ (conjugate, amber cap) into each well, mix gently.
3. **Cover the plate tightly and incubate for 3 hours at room temperature (18-24°C), in the dark.**
4. Aspirate and wash wells 5x with 300 µl diluted WASHBUF (wash buffer, natural cap). After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel.
5. Add 100 µl SUB (substrate, blue cap) into each well. Swirl gently.
6. **Incubate for 30 min at room temperature (18-24°C), in the dark.**
7. Add 50 µl STOP (stop solution, white cap) into each well. Swirl gently.
8. Measure absorbance immediately at 450 nm with reference 630 nm, if available.

## 8) CALCULATION OF RESULTS

Construct a standard curve from the absorbance read-outs of the standards using commercially available software capable of generating a four-parameter logistic (4-PL) fit. Alternatively, plot the standards' concentration on the x-axis against the mean absorbance for each standard on the y-axis and draw a best fit curve through the points on the graph. Curve fitting algorithms other than 4-PL have not been validated and will need to be evaluated by the user.

Obtain sample concentrations from the standard curve. If required, the pg/ml concentration can be converted into pmol/l by applying a conversion factor (1 pg/ml = 0.038 pmol/l; MW: 26 kDa). Respective dilution factors have to be considered when calculating the final concentration of the sample.

Example of a typical standard curve:

The quality control protocol supplied with the kit shows the results of the final release QC for each kit lot. Data for optical density obtained by customers may differ due to various influences and/or due to the normal decrease of signal intensity during shelf life. However, this does not affect validity of results as long as an OD of 1.50 or more is obtained for the standard with the highest concentration and the values of the CTRLs are in range (target ranges see labels).

### 9) ASSAY CHARACTERISTICS

<b>Method</b>	Sandwich ELISA, HRP/TMB, 12x8-well detachable strips		
<b>Sample type(s)*</b>	EDTA plasma, heparin plasma, citrate plasma		
<b>Sample volume</b>	50 µl / well		
<b>Assay time</b>	3 h / 30 min		
<b>Detection limit</b>	5.4 pg/ml		
<b>Standard range</b>	0 - 1600 pg/ml		
<b>Conversion factor</b>	1 pg/ml = 0.038 pmol/l; MW: 26 kDa		
<b>Precision</b>		<b>n</b>	<b>CV (%)</b>
	<b>Within-run</b>	3	≤8
	<b>In-between-run</b>		in progress
<b>Accuracy</b>		<b>n</b>	<b>Recovery (%)</b>
	<b>EDTA plasma</b>	6	100
	<b>Heparin plasma</b>	1	79
	<b>Citrate plasma</b>	1	103
	<b>Serum</b>	6	91

Dilution linearity of endogenous intact FGF23		n	Recovery of expected dilution (%)		
			1+1	1+3	1+7
	EDTA plasma	6	107	108	111
	Heparin plasma	1	99	97	107
	Citrate plasma	2	129	112	126
	Serum	6	87	74	67
<b>Specificity</b>	Endogenous and recombinant human intact FGF23				
<b>Use</b>	Research use only (RUO; in the process of CE registration)				
Values for apparently healthy donors		n	Median intact FGF23 [pg/ml]		
	EDTA plasma	22	24.4 pg/ml		
	Heparin plasma	22	26.4 pg/ml		
	Citrate plasma	22	17.4 pg/ml		
Serum	22	14.8 pg/ml			

\* This ELISA is optimized and validated for human plasma samples. Serum and other sample types are compatible with this ELISA. For further information on serum measurement and assay performance characteristics, matrix comparisons and stability data please visit our website [www.bmgrp.com](http://www.bmgrp.com) (see Validation Data) or contact our customer service by e-mail [info@bmgrp.com](mailto:info@bmgrp.com) or by phone +43/ 1/ 29107-45.

## 10) PRECISION

Within-run (intra-assay): 2 samples of known concentrations were tested 3 times within 1 kit lot by 1 operator.  
In-between-run (inter-assay): in progress

Within-run (n=3)	Sample 1	Sample 2
Mean (pg/ml)	99	800
SD (pg/ml)	8.4	10.2
CV (%)	8	1

## 11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date.
- Protect reagents from direct sunlight.
- Substrate solution should remain colorless until added to the plate.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

## 12) PRECAUTIONS

All test components of human source were tested against HIV-Ab, HCV-Ab and HBsAg and were found negative. Nevertheless, they should be handled and disposed as if they were infectious. All liquid reagents contain  $\leq 0.1\%$  Proclin 950 as preservative. Avoid contact with skin and mucous membrane. Proclin 950 is not toxic in concentrations used in this kit. It may cause allergic skin reactions – avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke, or apply cosmetics where reagents are used.
- Wear gloves, glasses, and lab coat while performing this assay.

- Sulfuric acid is irritating to the eyes and skin. Avoid contact with skin and mucous. Irritations are possible. Flush with water if contact occurs!

### 13) LITERATURE

1. Physiology of FGF23 and overview of genetic diseases associated with renal phosphate wasting. Bachetta J et al., *Metabolism* 2019; Epub ahead of print.
2. Renal and extrarenal effects of fibroblast growth factor 23. Vervloet M, *Nature Rev* 2019; 15:109-119.
3. Role of phosphate sensing in bone and mineral metabolism. Chande S and Bergwitz C, *Nat Rev Endocrinol* 2018; 14(11):637-655.
4. FGF23 Actions on Target Tissues-With and Without Klotho. Richter B and Faul CF, *Front Endocrinol (Lausanne)* 2018; 9:189.
5. FGF23 in Cardiovascular Disease: Innocent Bystander or Active Mediator? Stöhr R et al., *Front Endocrinol (Lausanne)*. 2018; 9:351.
6. FGF23 beyond Phosphotropic Hormone. Takashi Y and Fukumoto S, *Trends Endocrinol Metab*, 2018; 29(11):755-767.
7. Role of phosphate sensing in bone and mineral metabolism. Chande S and Bergwitz C, *Nat Rev Endocrinol* 2018; 14(11):637-655.
8. FGF23-Klotho signaling axis in the kidney. Erben RG and Andrukhova O, *Bone* 2017; 100:62-68.
9. Fibroblast Growth Factor 23-Mediated Bone Disease. Goncivlea AR and Jan De Beur SM, *Endocrinol Metab Clin North Am* 2017; 46(1):19-39.
10. Posttranslational processing of FGF23 in osteocytes during the osteoblast to osteocyte transition. Yamamoto H et al., *Bone* 2016; 84:120-130.
11. Coupling fibroblast growth factor 23 production and cleavage: iron deficiency, rickets, and kidney disease. Wolf M and White KE., *Curr Opin Nephrol Hypertens* 2014; 23(4):411-419.
12. FGF23 from bench to bedside. Kovesdy CP and Quarles LD, *Am J Physiol Renal Physiol* 2016; 310(11):F1168-1174.

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# # BI-20700 FGF23 INTACT ASSAY PROTOCOL AND CHECKLIST

## **PREPARATION OF REAGENTS:**

- Bring all reagents to room temperature (18-24°C).
- Prepare reagents and samples as instructed.
- Bring unused components to the storage temperature mentioned in the package insert.
- Take microtiter strips out of the aluminium bag and mark positions on the protocol sheet.

## **TEST PROCEDURE:**

- Step 1) Pipette 50 µl STD/CTRL/SAMPLE (standard/control/sample) into respective wells.
- Step 2) Add 50 µl CONJ (conjugate, amber cap) into each well. Swirl gently.
- Step 3) Cover tightly and incubate for 3 hours at room temperature (18-24°C), in the dark.**
- Step 4) Aspirate and wash wells with 300 µl WASHBUF (wash buffer, natural cap) five times. Remove remaining buffer by strongly tapping the plate against a paper towel.
- Step 5) Add 100 µl SUB (Substrate, blue cap) into each well. Swirl gently.
- Step 6) Incubate for 30 minutes at room temperature (18-24°C), in the dark.**
- Step 7) Add 50 µl STOP (stop solution, white cap) into each well. Swirl gently.
- Step 8) Measure absorbance immediately at 450 nm with reference 630 nm, if available.